

---

## Accelerate Detection With The 3DQuorum™, The Ultimate Breast Tomosynthesis Technology



---

Hologic is pleased to introduce [3DQuorum Imaging Technology](#) to the Breast Health Continuum of Care.

Designed to provide radiologists with fewer, thicker tomosynthesis slices to review, [3DQuorum](#) uses Genius AI™- powered analytics to reconstruct high-resolution 3D mammography data to produce 6mm SmartSlices, with 70-micron pixel resolution.

Each high-resolution SmartSlice overlaps the previous one by 3mm, ensuring there is no loss of 3D image data and continuity in scrolling.<sup>1-2</sup> Further to this, during the creation of SmartSlices, Genius AI identifies clinically relevant regions of interest and preserves important features.

### 3DQuorum represents our highest resolution 3D™ Imaging, faster

Compared to 1mm tomosynthesis slices, SmartSlices reduce the number of images to review by two thirds, giving an **average time-saving of one hour per day** in image interpretation time.<sup>1-2\*</sup>

#### References:

1. Hologic Data on File Technical report File: TFL-00059
2. Hologic Data on File Clinical Report: CSR-00116

Published on : Fri, 30 Oct 2020